




Laguba Healthcare

Heinz Oscar Gugger
Thiersteinallee 87
CH-4053 Basel


Telefon *41 (0) 61 331 20 50
Telefax *41 (0) 61 331 20 51
Natel *41 (0) 79 215 77 50

E-Mail: laguba@datacomm.ch
Homepage: www.laguba.ch
Bank: UBS Basel Kto.Nr. 11-125,014.0



Aqualizer™ Ultra

The Aqualizer™ Ultra is a new improved version of the Aqualizer™ that is designed for increased gum comfort and improved retention. Used for average size mouths. Our best seller.



Aqualizer™ Mini

The Aqualizer™ Mini is the new improved ultra shape with smaller pads and arch size. Used for kids and small mouths.



Aqualizer™ Classic

The Aqualizer™ Classic is our original design -- the first fluid-cushioned and balanced splint. The deep V in the front irritated some patients, so we designed the Ultra. As it is still preferred by some patients, we still make it available.

What makes the Aqualizer™ different?

The Aqualizer™'s revolutionary water system is different than any other product on the market. While most splints simply disable the bite long term and guess at optimal occlusion, the Aqualizer™ takes the guesswork out of treatment by allowing the body to naturally find functional balance.

The Aqualizer™ applies a physical law of nature called Pascal's Law, meaning that when you bite down on the Aqualizer™, the fluid is evenly distributed across the entire bite. The Aqualizer™ offers no opinion, enabling the body to naturally balance itself and function on its own.

Our bodies seek functional balance - with the Aqualizer™, it is found instantly:

The Aqualizer™ is the most accurate and effective differential diagnostic tool available. The Aqualizer™ allows the body to find instant functional balance by unraveling the distortion it is locked in, making it an invaluable tool to anyone who is suffering from the effects of [TMD](#). The Aqualizer™ obtains a perfect bite on the first try - splints made from this bite are immensely effective, predictable, and require remarkably few adjustments, saving valuable chair time.

Why the Aqualizer™?

The Aqualizer™'s unique water system:

- Immediately optimizes biomechanics
- Supports the jaw in comfortable position
- Removes the teeth from dominance, placing bite and body in harmony
- Straightens the bite to maximize other structures - Aqualizer™ is the bridge
- Normalizes the bite at will - without drugs
- Removes the imbalance trigger for patients who have existing malocclusion
- Enables systemic function and balance
- Allows the body to naturally balance itself
- Finds perfect occlusal balance, starting treatment immediately

Greater effectiveness for:

- Treating [TMJ/PDS symptoms](#), including headache, neck and shoulder pain and stiffness, and TMD related pain from trauma and whiplash. You can easily provide fast symptom relief at the first office visit.
- Orthodontics and cosmetics: The Aqualizer™ offers better functional treatment planning with a preview of harmony between teeth, mandibular position, TMJs, and the entire muscle complex. Relieves orthodontic triggered muscle pain during treatment, without inhibiting tooth movement.
- Oral surgery: Pre-surgical differential diagnoses and post-surgical pain and inflammation relief.
- Improves restricted opening problems when taking impressions: Wearing the Aqualizer™ prior to bite registration establishes a physiologic centric relation.
- Quick, accurate differential diagnosis: Use an Aqualizer™ as your first conservative, reversible, inexpensive test to definitively diagnose the problem.
- [Pain relief](#) with the Aqualizer™ confirms that the pain is of neuromuscular origin, triggered or exacerbated by the occlusion, rather than organic or intracapsular (joint) pain.
- Occlusal treatment planning: When occlusal treatment of any kind is indicated, have the patient wear an Aqualizer™ during the pre-treatment planning stage as a uniquely effective deprogrammer. This will allow the patient's own muscles to establish the mandibular position that produces maximum comfort and functionality. Incorporate this muscle compatible position into the treatment plan to get the most successful results.
- Build your splints to a muscle-guided occlusion: Step-by-step the Aqualizer™ precisely guides your adjustment of any occlusal surface, including rigid splints. Using an Aqualizer™ bite registration and/or deprogramming and tapping, the Aqualizer™ system pinpoints all interfering or displacing contacts so you can make that occlusal position permanent. Save adjustment time, eliminate trial and error, and increase patient comfort.

Rapid pain relief is the benefit most appreciated by more than one million patients. Health professionals safely and freely dispense the Aqualizer™ to any patient complaining of [head, neck, or shoulder pain](#). As a temporary painkilling treatment, the Aqualizer™ will effectively relieve neuromuscular pain that has its origin in the occlusion. [Contact us for more information on the Aqualizer™.](#)